

Self-retaining Vaginal Speculum

by Swolin

Made in polystyrene, single use



This speculum is opened and adjusted with one hand. It gives the operator a very effective exposure of cervix and the angle of the handle keeps the hand well away from the rectum. As it has only one lateral support for upper blade, the taking of smear samples etc. is easier. The two parts of the speculum are assembled and disassembled easily, but they cannot fall apart by themselves. The construction prevents the vaginal wall from being squeezed when the speculum is closed.

The Swolin speculum is delivered with a silicone band (latex free) around the blades during transport. This must be removed before use. It is delivered unsterile and is CE-marked.

It is packed 25 in a plastic bag. 8 bags are packed in a strong cardboard box containing 200 specula, which is the unit for shipping.

Ref. 12201